

Hybridization (overnight, @ 37 °C)

SAFETY: Hybridization buffer contains formalin, if spilled clean up and soak area in water, formalin extremely volatile and nasty, do not inhale

Cut $\frac{1}{8}$ of filter piece and place on slide according to which sample it is (3 filter pieces can go on a slide)

Slides are labeled with shapes and colors to help identify like so

Apply hybridization buffer and probe (9:1)

We multiply by two to get 18 μ L buffer and 2 μ L probe

Probe lives in small white box in the fridge labeled "ALV01", if you run out more aliquots

Can be found in -80 near atria. (Let Taylor or I know if there are less than 6 vials because we need to order it from Germany)

Buffer lives in -20 in lab on a shelf inside the door labeled "Hybridization Buffer"

If doing 6 samples: make 1 extra 18:2 buffer and probe

If doing 12 samples: make 2 extra 18:2 buffer and probe

Mix buffer and probe in a 1.5 mL Eppendorf tube

Apply 20 μ L to filter pieces along the "inner part of the pizza" - very viscous

Use the edge of a pipette to lightly spread the mix around to break the surface tension

Making hybridization chambers

Fold Kimwipe in fourths "hot dog" way and slide into 50 mL Falcon tube labeled on caps with "C"

Add slide with probe on top **make sure slide is horizontal at all times**

Add 1 mL of buffer in gap between slide and Kimwipe saturating the wipe

Place all the hybridization chambers in the incubation oven horizontally

To prevent from rolling around place a tray to keep tubes from rolling around

Wash buffer steps (2x @ 3 mL in wells @ 46 °C)

Mix Wash Buffer according to number of samples ***see Reagents for CARD-FISH***

Add 3mL per well

Gently remove the filter pieces handling on the crust of the filter and place into corresponding well

Place the saturated Kimwipe, slide, and uncapped hybridization chambers on Kimwipes in hood (formalin needs to evaporate and cannot be disposed of in trash until evaporated)

Place the wells with filters into the oven at 46 °C for 30 minutes

Transfer filters from first wash into second wash and replace in oven at same temp and time duration

TNT Buffer Wash (performed in wells for 15 minutes in dark at room temperature)

Add 3 mL of TNT Buffer into each well (check to make sure TNT Buffer is not cloudy or has precipitate, sign of bacterial growth)

Transfer filters to corresponding well and incubate for 15 minutes in the dark at room temperature

TSA and Fluoroscene (incubate in dark at room temp for 30 minutes)

Remove filters from TNT Wash and place on a slides (3 to a slide) in the previous order

Make sure the side of the filter with the sample is facing up!!

TSA and Fluoroscene lives in box in lab fridge

Add TSA and Fluoroscene system & amplification reagent in a 1.5 mL Eppendorf tube

For 6 samples: 2.5 μ L system & 125 μ L amplification reagent

For 12 samples: 5 μ L system and 125 μ L amplification reagent

Add 20 μ L per filter on slide in same arrangement as probe application and incubate in dark at room temp for 30 minutes

TSA and Fluoroscene system comes as a solid and needs to be dissolved in DMSO first, tell Taylor if running low and she will dissolve another vial. Keep an eye on reagent amounts so we can reorder

TNT Wash (2x at 55 °C)

Add 3 mL of TNT Buffer to each well

Place filters from TSA slides into the wells

Incubate for 20 minutes at 55 °C, repeat once in different wells with same amount of TNT Buffer

Calcafluor *stains cellulose in cells blue*

Remove the filters from the TNT wash and place on a slide (3 to a slide) with the sample side facing up

Add 20 μ L of Calcafluor to each filter piece

Calcafluor lives in a centrifuge tube in a beaker in a shelf in the lab fridge

Incubate in room temperature in the dark for 7 minutes

Millique Rinse (2x @ room temp for 1 minute each)

Fill wells with 3 mL Millique water

Remove the filters from the slides and place into wells, let sit at room temp for a minute

Transfer to next wash, wait for a minute

Citifluor/PI *Citifluor is a mounting reagent and PI (propidium iodide) stains the nucleus of the cell*

Remove filter pieces from the Millique rinse and place on final mounting slides with sample side up (**only 2 filter pieces to a slide**)

Add 20 μ L of Mounting Reagent (**refer to Reagents for CARD-FISH**) to each filter piece

Place coverslip over the filter making sure there is no air bubbles

Place mounted slides horizontally on a flat surface (well plates work well) in the lab fridge

Labeling slides

Label slides with corresponding shape and number, "Salt Pond", date of sample, cast and depth, filter pore size, volume filtered